



A Novel Biosensor for Water Quality Monitoring combines DNA Origami and Quantum Technology

Project GREENER

Christoph Meinecke (ENAS, TUC)



Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."

- Title: *"Single Photon source and detector based on novel materials for the detection of endocrine disruptors"*



GER: TUC (CO), ENAS, IFU, **SPA:** CICB, **FR:** AUREA, **IT:** SJ, WH, TC




- Health risks due to contamination (EDC) of water are a major problem
- Even extremely low concentrations of contaminants can have an impact on humans and animals
- Problem: Contaminants are not measurable on-site
- Target: On-site measuring of extremely low concentrations of hormones or hormone-like substances in water for fisheries and aquaponics
- Development of a **new type of single photon spectrometer**
 - Compact
 - Robust
 - Modular (source, detector, flow Cell with sample)




<https://albert-schweitzer-stiftung.de/wassertiere/fische-aquakultur>

<https://www.swd-ag.de/magazin/was-ist-aquaponik/>

- Focusing on:
 - 4,4'-Isopropylidendiphenol (Bisphenol A)
 - β -Estradiol
 - Estriol
 - Progesteron
 - 17 α -Ethynylestradiol
 - 3-Indolacetic acid (IAA)
 - Norethisteron
- NIR Reference measurements





β -Estradiol




Progesteron

Bisphenol A




NIR and MIR spectra,
substances dissolved in methanol @ 20mg/ml,
recorded with FTIR spectrometer VERTEX70

- Spectrometer and its Components
 - **Single Photon Source** (SPS): Array of SPS with different wavelengths → spectrum
 - **Flow cell**: Pre-concentration of target substances
 - **Single Photon Detector** (SPD): Improved single-photon detector with high efficiency in SWIR
 - **AI and Data cloud**: AI data analysis for improved selectivity
 - **Spectrometer set-up**: includes components, optical path adjustment, electronics



- Using **attenuated total reflection (ATR)** as measuring technique in conjunction with infrared spectroscopy
- Classical light sources (thermal, coherent)
 - High photon noise
 - Limits the low absorption measurements
- Using **single photon source** as light source combined with **single photon detector** with high quantum yield
 - High noise reduction
 - Measurement of substances in the ng/L range possible
- Development of Single Photon Source
 - Synthesis of new, green NIR-QDs integrated in DNA origami
 - Development and integration in LED-Layer stack



<https://www.leuchtmittelmarkt.com>



<https://germany.kyocera.com>


Brahim Lounis and Michel Orrit, *Rep. Prog. Phys.* **68**, 1129 (2005).

Ratio of variances:



$$\frac{\Delta T_{SP}^2}{\Delta T_C^2} = 1 - 2x \frac{T}{1 + T}$$

if $x \sim 1$ and $T \sim 1$:

→ perfect noise canceling



QD-LEDs combination with beam shaping elements




QD-LED

- Thin-film layer stack
- Components: Electrodes, HTL, ETL, QDs
- Injection of charge carriers causes radiative recombination

→ many photons emitted into half-space



Contacting of a single quantum dot using the QD-LED layer stack to realize a single photon source



Single Photon Source

- Contacting of one Quantum Dot
 - Using non-toxic, environmental friendly QD in SWIR range
- single photon emission

- Realizing precise localization of Qdot emitters
→ bottom-up DNA Origami Technology
- Folding of circular scaffold strand by shorter staple strands
- Staples overhang allows selective binding
→ Quantum Dots
- Due to the well defined position of each staples
→ nano bread board
- Resolution of 2,5 – 5 nm



• Immobilisation-Sequences (S3-S6)

S3: TCGATAATC-AGACTTTTTCATG

S4: ACGGTCGATA-GAGCGGGAGCTAAAC

S5: GTTTGACCATTAGA-AGTAGTAGAC

S6: ACGAGCTACC-TACATTTCGCAA

• NP-Functionalisation-Sequences (S1 & S2)

S1 (Top Left): CATGCATCAGGGTACTACGACGT-TGGGCGCATCAAACCC




S1 (Bottom Right): CATGCATCAGGGTACTACGACGT-TTACCATTAGCAA

S2 (Top Right): AAAGACAAAAGGG-CGATCGATCTCGAACGTACGACA




S2 (Bottom Left): TTTAACCAATAG-CGATCTGCTCGAACGTACGACA





- Fabrication of Qdot LED to realize the **Single Photon Source**
 - Wafer level technology for micro- and nanosystems
 - Deposition of functional layer to electrically contact the Qdot(s)
 - Preparing nanostructures to locally immobilize the DNA Origami
 - Deposition of DNA origami hybrids on nanostructured surfaces
 - Encapsulate the DNA-Qdot-hybrid
→ **(single) Qdot LED**



- Preparation of lenses on a transparent wafers
- Rear side preparation **(Single) Photon Source**
 - CdSe/ZnS QD-LED
- Measurement of the radiation pattern
 - Without lense: Lambertian radiator
 - With lense: narrow radiation




— 400 µm Wafer
— 500 µm Wafer
— without lense




Objectives

- Design of production of transparent fluid cell suitable for detection of analytes
- Development and optimization of methodology for preconcentration of analytes in vicinity to the optical detector
- Design of the analytical detection procedure



- Selective pre-concentration of analyte molecules for **label-free** detection



- Labeling analyte molecules with AuNP or enzymes

Andreas Morschhauser


andreas.morschhauser@zfm.tu-chemnitz.de

Martin Möbius

Project Coordinator
martin.moebius@zfm.tu-chemnitz.de

Julia Hann

Project Manager
julia.hann@zfm.tu-chemnitz.de



Susanne Hartmann

susanne.hartmann@zfm.tu-chemnitz.de

Mathis Janßen

mathis.janssen@zfm.tu-chemnitz.de

Jörn Langenickel

joern.langenickel@zfm.tu-chemnitz.de

Christoph Meinecke

meinecke@zfm.tu-chemnitz.de



Thanks for your attention



www.greenerproject.eu



Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."